



PRICE QUOTATION

CleanAir Express products are guaranteed to give you 100% satisfaction in every way. Return any CleanAir Manufactured product to us at any time if it proves otherwise. We will repair it, replace it, or refund the original purchase price. We do not want you to own any Clean Air Express product unless you are completely satisfied.

Date:

To:

Fax:

Regarding: CEM Heated 2 Stage Probe Quote

Quote No:

Thank you for your recent inquiry. We are pleased to submit the following quotation. Quotations are guaranteed for 30 days and are shipped F.O.B. from our facility in Palatine IL. Please contact me if you have any questions. Terms of payment are wire transfer or credit card for the first order. For order placement, fax the Purchase Order to (847)991-8924 or mail to the address above.

Part Number	Item	Price
07202	2' Effective CEM Heated Filter Probe with 20-micron SS Frit, Fiber Filter stem	\$721.00
07204	4' Effective CEM Heated Filter Probe with 20-micron SS Frit, Fiber Filter stem	\$784.00
07206	6' Effective CEM Heated Filter Probe with 20-micron SS Frit, Fiber Filter stem	\$832.00
07208	8' Effective CEM Heated Filter Probe with 20-micron SS Frit, Fiber Filter stem	\$890.00
07275	In-Line Self-Regulating (275 °F) Temperature Control (Option)	\$264.00

Our Field Emergency Service

CleanAir Express will ship, overnight and at our expense, any CleanAir Express component to complete your testing in the event of equipment failure including competing Method 5 equipment that is compatible with CleanAir components. We will do this for locations in the U.S. via FedEx Priority Overnight®. This service is provided free of charge to all Clean Air customers

Thank you,


Clean Air Engineering



CONFIDENTIALITY NOTICE

The documents accompanying this fax transmission contain confidential information. The information is intended only for the use of the individual(s) or entity named above. If you are not the intended recipient, you are notified that any disclosure, copying, distribution or the taking of any action in reliance on the contents of this faxed information is not permissible. If you have received this fax in error, please immediately notify Clean Air Engineering by telephone at the number above to arrange for the return of the original documents. If there are transmission problems or questions please contact Dan at the number above.

CEM Heated Filter Probe



The Clean Air Express CEM Heated Filter Probe has an integral filter chamber that will accept a glass fiber filter cartridge or a stainless steel sintered frit. Two stages of filtration may be achieved by using an in-stack coarse 20-micron SS frit on the probe end and a fine .02-micron filter in the chamber. The probe and filter chamber are heated by the same reliable high temperature heater eliminating the need for multiple temperature controllers. An optional in-line self-regulating (275 °F) temperature control thermostat may be incorporated into the power cord making the heated filter probe totally self contained and easy to use.

Voltage: 120/240 Vac

Part Number	Item
07202	2' Effective CEM Heated Filter Probe with 20-micron SS Frit, Fiber Filter stem
07204	4' Effective CEM Heated Filter Probe with 20-micron SS Frit, Fiber Filter stem
07206	6' Effective CEM Heated Filter Probe with 20-micron SS Frit, Fiber Filter stem
07208	8' Effective CEM Heated Filter Probe with 20-micron SS Frit, Fiber Filter stem
07275	In-Line Self-Regulating (275 °F) Temperature Control (option)